

Nothing Beats a Quickie®



rigid family



The Quickie® GP Series is the first family of rigid versions of the Quickie 2. The popularity of the GP family comes from the day-to-day versatility riders have come to rely on. The ideal chair for daily use and sport activities, the one-piece welded frame delivers a solid rigid ride.



improving people's lives

Quickie GP and GPV Manual Rigid Wheelchairs

Ultimate versatility for everyday and sports use



Innovative Design


Fully adjustable with a solid, rigid frame, Quickie GP and GPV chairs are engineered to maximize function and performance while minimizing effort for the rider. The extensive array of options and accessories make this model an ideal chair for daily use and sports activities. A variety of wheel options (tires, rims, hubs) are matched to the desires of the rider with choices from no maintenance to high performance, including Quickie Performance Wheels for tennis and basketball.



Family of Options

The GP family includes the Quickie GP, Quickie GPV, Quickie GPV Titanium, Quickie GPV Swing-Away and GPV Competition. The GP features a non-tapered front end, the GPV features a 70 degree V Front Frame and the GP Swing-Away has a rigid frame with swing-away footrests. The GP Titanium is the GPV made with titanium rather than aluminum and the GPV Competition is a package that offers an extra 2" of fixed tilt or "squeeze" for the rider.

With a variety of options, the GP series can be customized to fit tall, short, long or wide – any individual can find a comfortable home in a GP.



	GP (K0005)	GPV (K0005)	GP SA (K0005)	GPV Comp (K0005)
INNOVATIVE DESIGN				
Frame Style	Open-Traditional	Open-Traditional	Open-Traditional	Open-Traditional
Transport Weight	14 lbs	9.7 lbs	12.7 lbs	9.7 lbs
Standard Weight	26 lbs	21.5 lbs	24.5 lbs	21.5 lbs
Weight Capacity	250 lbs	250 lbs	250 lbs	250 lbs
CUSTOMIZATION/ADJUSTABILITY				
Seat Width	12" - 22"	12" - 22"	12" - 20"	14" - 20"
Seat Depth	12" - 22"	12" - 22"	12" - 22"	14" - 17"
Back Angle Adjustment	82° - 94°	82° - 94°	82° - 94°	82° - 94°
Camber	0° - 11°	0° - 11°	0° - 11°	0° - 11°
Tension Adjustable Uph	Yes	Yes	Yes	Yes
Overall Width	26.5"	26.5"	26.5"	26.5"
CLINICAL APPLICATIONS				
Front Seat-To-Floor Height	15.5" - 23"	15.5" - 23"	16.5" - 20.5"	17" - 21.5"
Rear Seat-To-Floor Height	14.5" - 20.5"	14.5" - 20.6"	16" - 20"	15" - 18"
Footrest/Frame Angle	60°, 70°, 75°, 80°	60°, 70°, 75°, 80°	Swing Away	80°
Back Type	Folding/Non Folding	Folding/Non Folding	Folding/Non Folding	Folding/Non Folding
Back Height	Adj 8" - 20"	Adj 8" - 20"	Adj 8" - 20"	Adj 11" - 20"
Center Of Gravity (COG)	.05" - 3.5"	.05" - 3.5"	.05" - 3.5"	.05" - 3.5"
LIFESTYLE OPTIONS				
Xtender-Power Assist	Yes	Yes	Yes	No
Front Suspension	No	No	No	No
Rear Suspension	No	No	No	No
Color Fusion	Yes	Yes	Yes	Yes
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime