



TRITON – FULL FRAME HEAVY USE DAY CHAIR



*The Triton is a handmade, lightweight, custom designed wheelchair build to the 'mm'. Available in chrome moly or titanium.
The Triton is designed for users who prefer a full frame wheelchair and also for those who are a little harder on their wheelchairs.*

STANDARD FEATURES

- Custom made chrome moly
- Lightweight
- ¾" diameter frame tubing
- Choice of frame colours
- Non folding backrest (fixed)
- Choice of adjustable or non-adjustable castor housing & rear wheel positions
- Adjustable height footrest
- Quick release spoke rear wheels C/with aluminium hand rims & high pressure or block tread tyres
- Choice of front castor wheel size
- Choice of scissor or push to lock brakes
- Velcro adjustable back & seat upholstery
- Integral push handles
- Calf strap

Transport Weight = 13 kgs

Frame Weight = 9.5 kgs

These weights are based on the following wheelchair configuration:

430mm X 420mm seat depth

Chrome moly frame

Fixed aluminium side guards

Aluminium adjustable footplate

24" Spoke quick release rear wheels comes with high pressure tyres & aluminium

Anodised push rims

4" Soft roll front castors