

The Snug Seat **GATOR** is a strong and stable gait trainer often chosen when a heavier more stable base is required. Available in a simple format of just the frame as a basic posterior gait trainer or with the addition of upper body support and leg separation as a full upright supported gait trainer. All accessories can be added or removed as skills change or improve.


The unique design of the Gator frame allows for an unrivalled amount of growth suitable for children aged 1½ to 7 years. This is accomplished with a height adjustable back frame, as well as the unique rotational handles, which also provide width adjustment.

Standard Features include: Easy-fold steel frame, height adjustable padded back bar, swivel front castors with pin locks, rotational hand grips and non-reverse brakes. Available in **Red** only.



Measurements

Width between handles	26.5 - 43cm
Height - floor to handles (K)	32 - 63cm
Width (G)	46cm
Length (L)	70.5cm
Weight	7.5 kg
Max. load/user weight	32 kg
Max. user height	115cm



Accessories



Side Supports
Provide hip and lateral support
Size Measurement
H: 11 x W:16 cm



Strap for back support
For added safety when walking
Size Measurement
31cm



Ankle Straps
Keeps child's legs and ankles in proper walking position
One Size



Strap for Side Supports
For added safety when walking
Size Measurement
L: 38-47cm



Back support on bar
Allows for maximum support of a child's trunk
Size Measurement
H: 10 x W:20 cm



Sling Seat
Provides partial weight bearing
One Size



Back Support
Allows for maximum support of child's trunk
Size Measurement
S H: 9 x W: 17 cm
L H: 10 x W: 24cm



Forearm support with hand grip
Perfect for upper extremity support and partial weight bearing
One size



Anti-Tips
Extra safety to prevent the Gator from tilting
One Size



Dejay Medical & Scientific Pty Ltd
1 Prince William Drive Seven Hills NSW 2147
Tel: +61 2 9838 8869 Fax: +61 2 9838 7869 Email: info@dejay.com.au
www.dejay.com.au