



Special Tomato jogger Push Chair - all terrain metro to park trails.

Perfect for on the go - a lightweight, maneuverable stroller that does not compromise the positioning.

With it's super-light aluminum frame, the Special Tomato™ Jogger folds with ease and can be placed neatly into your vehicle.

The comfortable padded seat reclines from 10°- 15° and there is a protective canopy with a window that allows you to keep an eye on your child.

For children who need minimal to moderate positioning, The Soft-Touch™ Liner or Sitter will easily attach.

Covered in easy-care, removable, washable, reversible contemporary fabric.

Standard Features

- Quick Release 12" pneumatic rear tyres
- Swivel lock 12" pneumatic tyre for easy maneuverability
- Foot well with side guard
- Functional rear suspension
- Height Adjustable Push handle with -
- attendant hand brake / front disc brake
- Teather Strap
- Rear parking brake-red
- Folding latch for quick and easy one step folding
- 5-Point harness
- Canopy with window - folds with jogger
- Removable Grab rail
- Reclining Seat Back
- Latex-Free padded seat and back surface



Specifications

Seat to Foot Well:

Stroller Only	35 cm
Liners Size 1 & 2	40 cm
Sitter Size 1	42.5 cm
Sitter Size 2	45 cm

Seat:	35w x 31d cm
Back Height	60 cm
Max. User Weight	50 kgs
Stroller Weight	12.3 kgs

Dimensions

Overall	66w x 107l x 110h cm
Folded (without rear wheels)	54 - 84 cm
Handle Height to Floor	89 - 110 cm
Seat can be set at	50°



Pictured Above:

Jogger shown with optional Soft Touch Sitter **Chocolate** with **Lilac** Harness

Jogger shown with optional Soft Touch Liners in **Cherry**

Optional Accessories

- Padded Head Support
- Medical Necessities Bag

Size one Soft Touch Liner Seat and/or Back Cushion

Size two Soft Touch Liner Seat and/or Back Cushion

Soft Touch Liners available in **Cherry**, **Chocolate** and **Lilac**

Size one Soft Touch Sitter with Harness - available in **Chocolate** only

Size two Soft Touch Sitter with Harness - available in **Chocolate** only

