



**SPECIAL
TOMATO**

Soft-Touch Liners

Special Tomato Soft-Touch Liners - Stylish, Durable, Portable, Seating

Adding comfort to your child's sitting experience is very important!

Special Tomato has used its unmatched technology to create the Soft-Touch Liners. These liners offer support and comfort for your child with mild to moderate physical involvement.

Soft-Touch Liners Size 1 and 2 can be easily inserted into the Special Tomato EIO Push Chair and or Jogger (pictured below). All three sizes can be attached to any standard chair found in the home, school or day care.



Available in three delicious colours **Cherry**, **Chocolate** and **Lilac**

The liners are made from Latex-Free, antimicrobial soft-to-the-touch materials that won't tare or puncture, impermeable to fluids and can be safely wiped clean with mild disinfectants - perfect for use with multiple children.

Special Tomato Soft-Touch Liners - are the accessory you have been looking for to ensure comfort for your child.

Standard Features

- Latex Free
- Anti-microbial built in protection
- Impermeable to fluids
- Peel and Tear Resistant
- Contoured Seat and Back
- Built in Abductor
- Attachment Straps
- Can be used together or separately
- Available in 1,2,3 different sizes and



Pictured Above: Soft-Touch Liners in optional Special Tomato EIO and Jogger

Special Tomato Soft-Touch Liners Sizing Chart

	Size 1	Size 2	Size 3
Seat Cushion Width	20 cm	25 cm	30 cm
Seat Cushion Depth	20 cm	25 cm	30 cm
Overall Dimensions of Seat Cushion	27d x 28w cm	33d x 31w x 10h cm	38d x 36w x 11.25h cm
Overall Height of Seat Cushion	8.75h cm	10h cm	11.25h cm
Weight of Seat Cushion	.68 kg	.91 kg	1.13 kg
Back Cushion Width	18.75 cm	21.25 cm	25 cm



Sydney Head Office - Ph: (02) 9838 8869 Email: info@dejay.com.au
Victorian Office - Ph: (03) 9558 9669 Email: jackie@dejay.com.au - Queensland Office - Ph: (07) 3863 2589 Email: malcolm@dejay.com.au

Importers & Manufacturers of Rehabilitation Mobility Equipment