



**SPECIAL  
TOMATO**

# Soft-Touch Sitter


**Special Tomato Soft-Touch Sitter**– soft, supportive, comfortable seating experience for children with special needs.

Available in 1,2,3,4,5 different sizes and 3 different systems, the Soft-Touch Sitter is the ideal choice for adaptive seating.

This soft-to-the-touch seat provides a seating experience for a special needs child that is comfortable and functional.

The Soft-Touch Sitters are made from a Latex-Free, seamless material that is impermeable to fluids and guaranteed to be peel and tear resistant. Special Tomato Soft-Touch Sitters have a strong exterior that will stand up and last, while the inside is a soft, comfortable foam, size's 1 and 2 can be easily inserted to the Special Tomato EIO or Jogger.

Available in Chocolate only with a choice of **Lilac** or **Tan** harness



## Soft Touch Sitter Standard Features

- Contoured head and lateral supports
- Soft contoured seat with anti-thrust seating surface
- Built-in abductor
- Adjustable 5-Point positioning harness
- Attachment straps for securing to a Tomato base or standard chair

## Choose from three Systems

- Soft-Touch Sitter Only
- Soft-Touch Sitter with Wooden Tilt Wedge Kit
- Soft-Touch Sitter with Wooden Mobile Tilt Wedge Kit

## Soft-Touch Sitter with Wooden Tilt Wedge Kit - Wooden Mobile Tilt Wedge Kit

Offering proper seating positioning for children with mild to moderate physical involvement, the Soft-Touch Sitter with Tilt Wedge Kit lifts the Soft-Touch Sitter off the ground and brings children with special needs to the level of their peers while playing on the floor.

With a 25° of tilt the Special Tomato Soft-Touch Sitter can be positioned on the Tilt Wedge in several different angles, providing the perfect support for your child.

With the optional Wooden Mobile Tilt Wedge Kit, includes four swivel casters with rear locking brakes, allowing your child to move from one area to another with ease and without compromising their posture and positioning.



## Special Tomato Soft-Touch Base Sizing Chart

	Small		Large	
	Stationary Wood Tilt Wedge	Mobile Wood Tilt Wedge	Stationary Wood Tilt Wedge	Mobile Wood Tilt Wedge
Use With Sitter Sizes	1, 2, 3	1, 2, 3	4, 5	4, 5
Over width	17.5cm	17.5cm	20.5cm	20.5cm
Overall Depth	30cm	30cm	35cm	35cm
Overall Height	16cm	19.5cm	23cm	26.5cm
Weight of Base	12kg	14kg	15kg	17kg



## Special Tomato Soft-Touch Sitter Sizing Chart

	Size 1	Size 2	Size 3	Size 4	Size 5
Target Weight	9-14 kg	14-27 kg	18-36 kg	23-50 kg	41-91 kg
Target Height	75-87.5 cm	87.5-120 cm	100-140 cm	125-162.5 cm	135-185 cm
Target Age	1 to 3	3 to 6	6 to 9	9 to 14	13 to Adult
Inner Seat Depth	14.5 cm	23.75 cm	28.75 cm	33.75 cm	40 cm
Inside Seat Width		25 cm	30 cm	35 cm	40 cm
Seat Back Height	52.5 cm	65 cm	72.5 cm	82.5 cm	97.5 cm
Overall Seat Depth	30 cm	35 cm	40 cm	47.5 cm	52.5 cm
Overall Seat Width	30 cm	35 cm	40 cm	45 cm	50 cm
Overall Seat Height	62.5 cm	72.5 cm	85 cm	95.25 cm	110 cm
Weight of Seat	2.7 kg	4.5 kg	5.9 kg	7.25 kg	8.6 kg
Floor to Seat Ht UPRIGHT on Stationary Base	23.5 cm	23.5 cm	24 cm	23-50 kg	41-91 kg
Floor to seat Ht TILTED on Stationary Base	23.5 cm	24.5 cm	25 cm	33 cm	36.2 cm
Floor to seat Ht UPRIGHT on Mobile Base	32 cm	33 cm	33 cm	37.5 cm	40.6 cm
Floor to seat Ht TILTED on Stationary Base	32.5 cm	34.5 cm	34.5 cm	40.6 cm	45.1 cm



Sydney Head Office - Ph: (02) 9838 8869 Email: [info@dejay.com.au](mailto:info@dejay.com.au)  
 Victorian Office - Ph: (03) 9558 9669 Email: [jackie@dejay.com.au](mailto:jackie@dejay.com.au) - Queensland Office - Ph: (07) 3863 2589 Email: [malcolm@dejay.com.au](mailto:malcolm@dejay.com.au)

Importers & Manufacturers of Rehabilitation Mobility Equipment