

Find Your Freedom



With Our Adaptive Bicycles

Adjustable Features That Get *Everyone* Riding

Non-intrusive Rear Steering™

- Allows caregiver to steer and brake from back of bike
- Optional jogging extension handle with on-the-fly neutral gearing



Adjustable Handlebars

- Quick release lever allows handlebars to move out of the way for transferring
- Fully-adjustable handlebar height and angle



Seating and Headrests

- Custom to any ability
- Height-adjustable
- Best support in the industry
- Multiple headrest options available



Lock/Unlock Steering

- Standard on all models
- Provides the option of rider steering or caregiver steering



Adjustable Tilt Fork Accessory

- Add or reduce the amount of seat tilt for perfect positioning
- 6" of tilt range, translating to a variable 10° tilt



Custom-Fitted Foot Plates

- Variety of foot plate options
- Available with abductor leg braces
- Assists with leg and knee alignment



Tires and Rims

- Available with solid or pneumatic tires
- Optional heavy-duty rims



**Narrow and Split Frames Available on Select Models for Indoor-Outdoor Use and Easy Transportation*


Customize Your Bike

All Freedom Concepts adaptive bicycles are highly adjustable and built to last. If none of our standard bikes work for you, we can custom build something that will. Contact us today to learn more.

freedomconcepts.com • 1-800-661-9915

Electric Assist

- Excellent for longer rides
- Assists with uphill pedaling
- Long lasting battery life
- Teaches independent riding



Lithium-Ion Battery



Hub Motor

"The electric assist motor has made all of the difference. It allows us to go for longer rides as a family, and Ben absolutely loves it. He weighs about 110 pounds, so when you add the weight of the bike, I'm not strong enough to push him up hills. With the electric assist, Ben can pedal up hills all by himself. We wouldn't be able to ride it without it."

- Courtney Patton, Ben Patton's mother.



Adjustable handlebars leave lots of space for transfers.



Height adjustable seating systems make for smooth transfers everytime.

Transfer With Ease



Made To Last

- Highly Adjustable To Accomodate Growth
- Custom Built In-House
- Thoroughly Inspected Before Shipment



For more information, visit freedomconcepts.com or talk to your Dejay Medical Representative

Follow Us!    @freedomconcepts